

MACH LED 130 CEILING MOUNTED LIGHT WITH 40CM CEILING KIT - CLEARANCE

MACHLED130-C/KIT-40-CL



PRODUCT DESCRIPTION

MACH LED 130 Focusable Ceiling Light - With Lamp Head - 40cm

Mach LED 130 Focusable Light & Ceiling Kit - Pole length set at 40cm

Consists of:

- Ceiling Tube
- Ceiling Flange
- Swing Arm
- Mach LED 130 Lamp Head For Ceiling Mounted Light

Approx ceiling height of 2.8m

*Please note: Anchorage Ring is not included

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED130-C/KIT-40-cl
Brand	Dr Mach

Autoclavable Handle

Yes

Focusable

No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights