

## MACH LED 130 EXAMINATION CEILING LIGHT

MACHLED130-C-CONFIG\*



### PRODUCT DESCRIPTION

## The Mach LED 130 is a compact examination light with powerful 60,000\* lux LED technology

Producing a cool light the LED technology is much more effective than conventional light sources. The temperature increase in the surgeons working area is almost non-existent. Features include a faceted multi lens system, a colour rendering index of Ra=95 and superior handling design allowing the light to be positioned exactly to the wound field.

Available with optional focusable head allowing easy viewing of the working field and sterilizable handle

When the ceiling height is below 2.8m or above 3.0m please state the fixation height and height to underside of the false ceiling so the correct length pole can be supplied.

To calculate pole length required, deduct 2400mm from your ceiling height.

#### Key Features:

- 60,000 Lux (at a distance of 1m)
- Faceted multi lens system
- Colour rendering index of Ra=95
- Authentic colour rendering to reduce eye strain
- Superior handling design - allowing light to be positioned where needed
- 0.5°C temperature increase
- Focusable
- Colour temperature (Kelvin): 4300

- Dimming range (%): 50-100
- Number of LED: 19
- Size of the light field (cm): 12 (fix)
- Total power consumption (W): 28
- Life-span of the LED's (h): 50.000
- Working distance (in cm): 70-140
- Chassis of the light (cm): 33

## Hygienic

High quality materials provide security and avoid infiltrations of dust, dirt and liquids inside the lamp head. The disk sealings of the light outlets and the circumferential sealing cord guarantee a perfect product down to the smallest detail.

## Effective

Powerful LED-technology helps keeping a cool head through an almost non-existent heat development.

## ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	MACHLED130-C-Config*
<b>Colour/Style</b>	White
<b>Brand</b>	Dr Mach
<b>Battery Backup</b>	No
<b>Colour Control</b>	No
<b>Colour Temperature (k)</b>	4300
<b>Head Type</b>	Single
<b>LED Life (hours)</b>	50,000
<b>Light Source</b>	LED
<b>Lux @1 meter</b>	60000
<b>Produceure</b>	Dental, Exam, Minor-Operation

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

