

## BURTONS VETSPIN DUO CENTRIFUGE

147-300



### PRODUCT DESCRIPTION

**The VetSpin Duo has been designed to meet the specific needs of a busy veterinary practice, with 3 quick pre-set programs for ease and speed of use.**

Its combination rotor accepts capillary tubes, microtubes and blood tubes making it ideal for simultaneously processing blood samples, semen, urine, DNA & RNA, organelles, microbial cells, and more.

An ergonomic design allows for rapid testing.

#### Key Features

- Range of 3 pre-set programs with the option to adjust the time and speed.
- Reaches high speeds of RCF 9,500 with speed accuracy of 100rpm for quick and accurate use.
- Solid and durable quality; fitted with suction stands for additional stability.
- Combination rotor fit for varying microtubes.
- Low noise characteristics of below 65dB.
- 1-99-minute timer setting.
- Components that will withstand years of hard service.
- Lid lock safety system.
- Imbalance switch prevents damage to centrifuge.
- 2 open lid detectors: 1 magnetic and the other on the lid-lock. The centrifuge will not run with an open or compromised lid.

#### Specifications

- Rotor: Combination Rotor for Haematocrit tubes, microtubes, blood tubes & 75ml capillary tubes
- Speed: 10,000rpm
- Dimensions: 208mmH x 305mmL x 270mmW
- Mains Frequency: 50/60Hz
- Fuses: 3amp 230V fuse
- Motor: AC Brushless Motor
- Timer: 0-99-minute timer
- Brake: Automatic
- Mains voltage: 220-230V
- Speed accuracy: 50 rpm
- RCF: 9,500g max
- Weight: 7.5kg
- Noise Levels: Below 65dB

**The box will include:**

- 1 x VetSpin Duo Combination Centrifuge
- 1 x Combination Rotor (for capillary tubes, microtubes, blood tubes)
- 1 x User Manual
- 1 x Rotor Removal Tool Set
- 1 x Lid Cushion
- 1 x 3amp Fuse
- 1 x Rotor Lid
- Warranty Card
- Starter kit including 2 x vials of capillary tubes and clay 1 sheet of clay

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	147-300

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights