

WHITE ECG ELECTRODE PADS - 3.6CM WITH POPPER X 50

VSL-WS00S50

**PRODUCT DESCRIPTION**

The Ambu WhiteSensor WS features a highly conductive solid gel with strong adhesion to ensure a good signal quality during short-term ECG monitoring applications.

Thanks to the flexible foam backing material, the electrode is easy to apply.

Key Features:

- Highly conductive wet gel - immediate reliable recordings
- Comfortable foam backing - resistant to fluids
- 3 or 4 mm fitting connector - direct connection to cable, without adapters
- Radiolucent and MRI Conditional - can remain on the patient during X-ray and MRI procedures
- Convenient presentation of 5 pieces/strip - for fast application

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	VSL-WS00S50
Monitor Fit	Vetronic Lightning, Vetronic Sapphire Plus, Vetronic Sentinel
Accessory Parameter	ECG

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights