

ECG ELECTRODE PADS, 3.5CM ROUND, WITH POPPER (X50 PACK)

VSL-WS00S50



PRODUCT DESCRIPTION

The Ambu WhiteSensor WS features a highly conductive solid gel with strong adhesion to ensure a good signal quality during short-term ECG monitoring applications.

Thanks to the flexible foam backing material, the electrode is easy to apply.

Key Features:

- Highly conductive wet gel - immediate reliable recordings
- Comfortable foam backing - resistant to fluids
- 3 or 4 mm fitting connector - direct connection to cable, without adapters
- Radiolucent and MRI Conditional - can remain on the patient during X-ray and MRI procedures
- Convenient presentation of 5 pieces/strip - for fast application

Change as used. Compatible with Edan X12, Mindray uMEC, Mindray ePM12, BLT E30VET & SurgiVet Advisor monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	VSL-WS00S50
Brand	Vetronic
Accessory Parameter	ECG
Monitor Compatibility	BLT E30VET, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor
Change Frequency	As Used

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights