

## 2.3MM OVAL SNARE - 30MM X 2350MM L (END-1100) - CLEARANCE

END-1100-CL



### PRODUCT DESCRIPTION

---

#### Flexible Oval Snare 30mm - 2.3mm D x 2350mm L

Highly reliable and easy to maintain, our flexible endoscopic instruments offer consistent performance and are CE and FDA approved.

- Constructed of special corrosion resistant steel alloys
- Autoclavable handles
- Colour coded and marked for identification

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	END-1100-CL
<b>Procedure Type</b>	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
<b>Channel Size</b>	2.3mm
<b>Working Length</b>	2350mm

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights