

## NON-AUTOCLAVABLE WATER BOTTLE FOR OLYMPUS ENDOSCOPES - CLEARANCE

END-WB-100-CL



### PRODUCT DESCRIPTION

---

## Generic replacement water for the 100, 130 Series Olympus Endoscopes.

This replacement water bottle has a 250ml capacity and offers the following features:

- Easy to connect/disconnect from scope
- Metal connector to bottle top and metal endoscope connector for durability
- Attached cap to protect endoscope connector
- Metal clip and strong plastic strap to attach the water bottle to the stock hanger on Olympus light sources
- Non-autoclavable - use soaking methods to sterilise

### ADDITIONAL INFORMATION

**Animal**

N/A

**SKU**

END-WB-100-cl

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights