

MOTIC 2820 LED MICROSCOPE - CLEARANCE

145-350-CL



PRODUCT DESCRIPTION

The Motic 28 series offers a unique ergonomic shape as well as superior optics than the standard DIN objectives.

A reliable, high-quality microscope at an affordable price.

Key Features

- Binocular head, inclined 45° and 360° rotating
- Widefield eyepieces WF10X / 18mm
- Quadruple nosepiece
- Achromatic objectives 4x, 10x, 40x, 100x (oil)
- Coaxial coarse and fine focusing adjustment
- Built-in low positional coaxial mechanical stage
- Rack & pinion focusable 1.25 N.A. Abbe condenser
- Iris diaphragm with filter holder
- LED illumination with intensity control for truer tissue colour compared to traditional halogen

Technical Specifications

- Optical System: Finite optical system, 160mm
- Observation tube: Monocular
- Interpupillary Distance: 55-75mm
- Diopter adjustment: On both tubes, +/- 5 diopter
- Inclination: 45 inclined, 360 rotating
- Eyepiece(s): WF 10X / 18mm

- Nosepiece: Reversed quadruple revolving nosepiece with click stops
- Objective classification:EA Achromatic, DIN
- Objective information: 4X / 0.10, 10X / 0.25, 40X / 0.65 / S
- Objective mounting thread: W 4/5" x 1/36" (RMS standard)
- Stage: Built-in low position coaxial mechanical stage with sample holder
- Stage size: 133 X 136mm
- Mechanical stage X & Y range: 77 x 31mm
- Upper limit stop: Preset but adjustable
- Condenser:Rack & pinion focusable 1.25 N.A. Abbe Condenser with iris diaphragm
- Focus mechanism: Coarse and fine focusing system with tension adjustment
- Fine focus precision: 2um minimum increment
- Z-axis movement: 12mm
- Filter: Built-in frosted glass
- Illumination: LED 20mA, 3.5V, 70mW with intensity control
- Transformer: External
- Power supply: Charger 100 - 240V or rechargeable batteries
- Accessories included: Dust cover, charger and batteries

S/N: 03240002784

ADDITIONAL INFORMATION

Animal	N/A
SKU	145-350-cl
Overall Product Dimensions (mm)	215L x 172W x 380H
Brand	Motic
Lens System	Binocular
Bulb Type	LED
Objective Type	Achromatic
Digital Camera	No
Supplied Objectives	100x (Oil), 10x, 40x, 4x
Microscope Complete Kit	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

