

ELITE RETINOSCOPE STREAK

WA18245



PRODUCT DESCRIPTION

The Welch Allyn Elite Retinoscope makes refractions easier than ever by producing one of the brightest, sharpest, and cleanest streaks available.

Enhanced light output makes refractions easier than ever, by producing one of the sharpest, brightest, and cleanest streaks available on the market. Discover unparalleled performance in a compact design with convenient, user friendly features.

Key Features

- Provides 2.8X the brightness of our standard retinoscopes
- Superior streak quality enables easier and quicker refractions even in patients with small pupils
- Crossed linear polarizing filter eliminates bothersome glare from trial lenses, helping you perform a more precise exam
- External focusing sleeve improves comfort and maneuverability
- Sealed optics keep the lenses cleaner for more effective exams
- Fits all Welch Allyn 3.5 V power sources (Handle Not Included)

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	WA18245
Overall Product Dimensions (mm)	17.2cmW x 4.7cmH x 4.7cmD
Generic Size	X-Small

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights