

## VACUUM VETCLAVE - CLEARANCE

2545-B S3 DW A-CL-3



## PRODUCT DESCRIPTION

**Designed specifically for the needs of a busy veterinary practice. The Burton's Instaclave range has been an industry favourite for over 30 years.**

The Instaclave 2545 Series 3 Type B is a Vacuum Autoclave designed specifically for the decontamination of wrapped porous loads and surgical instruments. It still has the added benefit of offering a quick non-vacuum cycle option, making it ideal for all instrument sterilisation in your busy surgery.

With a front-loading chamber, a water replenishment process that's easily accessible and cycle selection at the touch of a button, it is compact and designed to last but most importantly it is user-friendly and efficient.

Choose from a non-vacuum rapid cycle, a vacuum cycle for solid instruments or a cycle specifically for wrapped porous loads. This model has two sterilisation temperature options in both non-vacuum and vacuum modes. These can all be selected by the push of a button, with no complicated touchscreens, and therefore, reducing the number of components that can go wrong and cause maintenance downtime.

Unlike some other larger vacuum autoclaves in the industry, this autoclave is 600 mm wide which allows this machine to be placed on a benchtop within the practice, taking up a limited workspace, but allowing for adequate airflow. The large chamber size of 22-litre capacity enables more instruments to be sterilised at any one time. The Instaclave has space for 6 trays. 3 trays are included as well as your hose, drip tray and handle for ease of use. The door gasket can be replaced and the overheat cut out can be reset easily by the user, reducing your maintenance costs and downtime.

### What Instruments Can It Sterilise?

- Solid Surgical Instruments Wrapped
- Solid Surgical Instruments Unwrapped
- Wrapped Porous Loads

## Integrated Safety Features

You will always feel comfortable using the Instaclave due to its abundance of safety features. The manual and mechanical door locking mechanism ensures that it cannot be opened whilst sterilising due to the pressure forcing the second lock to engage. Overheating protection and a pressure safety valve are two other features that will guarantee you an easy and safe experience.

To be compliant with European Standards, Vacuum sterilisation really is the best way forward when choosing your new autoclave. [Read our guide on why a vacuum autoclave is an essential piece of equipment in some practices.](#)

## Key Features:

- 3 Fractionated Air Removal Stages
- 4 Cycle selection types
- Large chamber and tray size
- Safety features including manual and mechanical door locks and overheating protection
- Easy to use, clear buttons to switch between cycles, with clear lights to show your selection
- Non Drying cycle option
- Fully Integrated clean and dirty water reservoirs
- Software that allows the user-visible cycle pass or fail indication
- User Friendly Daily/Weekly maintenance checks

| Technical Specifications |  |
|--------------------------|--|
| Cycle Types              | 134°C Full Fractionated Vacuum Cycle<br>121°C Full Fractionated Vacuum Cycle<br>134°C Full Fractionated Vacuum No Drying Cycle<br>121°C Full Fractionated Non-Vacuum Cycle |
| Chamber Size             | 22 Litres  |
| Capacity                 | 6 Tray (Supplied With 3)   |
| Tray Size                | 400 mm x 180 mm x 18 mm DWH  |
| Dimensions               | 540mm x 500mm x 410mm WDH  |
| Cycle Times              | 134°C Full Fractionated Vacuum Cycle (80 mins)<br>121°C Full Fractionated Cycle (80 mins)<br>134° Full Fractionated Vacuum No Drying Cycle (30 mins)                       |

**S/N: 055484**

## ADDITIONAL INFORMATION

|  |                         |
|--|-------------------------|
| <b>SKU</b>                             | 2545-B S3 DW A-cl-3     |
| <b>Overall Product Dimensions (mm)</b> | 540mm x 500mm x 410 WDH |

|                                   |                 |
|-----------------------------------|-----------------|
| <b>Generic Size</b>               | Medium          |
| <b>Colour/Style</b>               | Stainless Steel |
| <b>Number of Cycles</b>           | 4               |
| <b>Vacuum</b>                     | Yes             |
| <b>Chamber Size</b>               | 22 Litres       |
| <b>Autoclave Accessory/Spares</b> | Yes             |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights