

TERUMO NEEDLE BOX OF 100 - CLEARANCE

400-0-NEEDLES-CL

**PRODUCT DESCRIPTION**

Terumo's needles use a unique double bevel design and precision manufacturing to provide ultra sharp needles, which translate into greater comfort for the patient

This product contains no components made of natural rubber latex and their tamper-evident blister packs result in the minimal waste disposal

- Terumo Needle, supplied in a box of 100
- UltraSharp Needles
- Tamper Evident Peel Packs

ADDITIONAL INFORMATION**Animal**

N/A

SKU

400-0-needles-cl

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights