

X-RAY & ULTRASOUND FOAM POSITIONER

END-USXP1



PRODUCT DESCRIPTION

Ensure a comfortable yet secure procedure for your cats and dogs with this alternative positioning aid. Suitable for most body weights and sizes in everyday general practice

- Dimensions: 100 cm x 35 cm x 14 cm + 2 x 70 cm long holding straps
- Weight: 2.3 kg

"I have found this product irreplaceable. This level of patient comfort and positional security is unobtainable using rigid plastic or short foam restraining aids - it is unique. Even with light sedation, patients will lie in this trough in comfort due to feeling secure and padded. This enables a faster and more accurate imaging procedure."

Dr Sophie Betts

Northeast Veterinary Imaging Ltd

ADDITIONAL INFORMATION

SKU

END-USXP1

Animal

Small Animal

Overall Product Dimensions (mm)

1000 x 350 x 140

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights