

VETRONIC SAV AIR PUMP

VSL-VS0460



PRODUCT DESCRIPTION

Want to ventilate your reptiles on room air but find it too fiddly to set up? The SAV Air Pump is a quick, affordable and easy solution.

- Affordable
- Efficient
- Perfect for your Reptilian and Avian patients
- 3 Pre-set Flow Rates

The SAV Air Pump is a reliable means of supplying room air to your SAV machine where ventilation with 21% oxygen is required for long or short- term support.

Preset flow rates and/or a controlled leak allow the use of this pump in all size animals up to 10kg.

Supplied ready to connect & use!

ADDITIONAL INFORMATION

| | |
|---------------|------------|
| Animal | Mixed |
| SKU | VSL-VS0460 |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights