

**EDAN HANDHELD PULSE OXIMETER CHARGING BASE**

E-VEH100B-POD

**PRODUCT DESCRIPTION****Edan Pulse Oximeter Battery Charger Stand Kit (UK)**

Kit Includes:

- 1 x Ni-MH Battery
- 1 x Charger Stand with Power Cord

**1 year warranty****ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	E-VEH100B-POD
<b>Monitor Parameters</b>	SpO <sub>2</sub>
<b>Screen Type</b>	Digital Display
<b>Power Source</b>	Battery, Rechargeable

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights