

## EDAN HAND HELD PULSE OXIMETER CHARGING BASE

E-VEH100B-POD



### PRODUCT DESCRIPTION

Edan Hand Held Pulse Oximeter Charging Base. Change when broken. Compatible with Edan VE-H100B.

## Edan Pulse Oximeter Battery Charger Stand Kit (UK)

Kit Includes:

- 1 x Ni-MH Battery
- 1 x Charger Stand with Power Cord

**1 year warranty**

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	E-VEH100B-POD
<b>Accessory Parameter</b>	Other
<b>Monitor Compatibility</b>	Edan VE-H100B
<b>Change Frequency</b>	When Broken

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights