

HEINIGER BLADE, SIZE 50

HNG-707-905-A

**PRODUCT DESCRIPTION****Heiniger Snap-On #50 / 0,2 mm Blade**

For veterinary use only. Interchangeable blade set for Heiniger Opal.

Heiniger clipper heads are compatible with clippers from Andis, Moser, Oster and Wahl.

0.2mm

ADDITIONAL INFORMATION

| | |
|---|---------------|
| Animal | Small Animal |
| SKU | HNG-707-905-A |
| Clipper Speed(s) | N/A |
| Number of Clipper Speeds | 1 |
| Compatible with other brand blades ? | Yes |
| Replaceable Battery | N/A |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights