

**HEINIGER BLADE, SIZE 40**

HNG-707-910-A

**PRODUCT DESCRIPTION****Heiniger Snap-On #40 / 0.25 mm Blade**

For veterinary use only. Interchangeable blade set for Heiniger Opal.

Heiniger clipper heads are compatible with clippers from Andis, Moser, Oster and Wahl.

0.25mm

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	HNG-707-910-A
<b>Clipper Speed(s)</b>	N/A
<b>Number of Clipper Speeds</b>	1
<b>Compatible with other brand blades ?</b>	Yes
<b>Replaceable Battery</b>	N/A
<b>Blade Size</b>	40

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights