

PANOPTIC OPHTHALMOSCOPE SET INCLUDING LITHIUM ION HANDLE CASE & CHARGER - CLEARANCE

WA11824-VSM-CL



PRODUCT DESCRIPTION

Thanks to an innovative, proprietary optical design known as Axial Point Source Optics, this ophthalmoscope allows all veterinarians to easily examine the fundus.

- Fast, easy entry into small, undilated pupils
- 5X larger view of the fundus vs. standard ophthalmoscopes in an undilated eye
- 25° field of view versus the 5° field of view of standard ophthalmoscopes
- Greater working distance improves comfort for practitioner
- Six inches from veterinarian to patient's eye
- Complete compatibility with existing Welch Allyn 3.5 V power sources

Patient Eyecup: The eyecup establishes and maintains proper viewing distance, provides image stability, serves as a pivoting point for panning around the retina, helps with orientation, and screens out ambient light. The eyecup is cleanable and reusable.

PanOptic™ Soft Grip: Ergonomically designed for maximum comfort, balance, and access to controls.

Patented Axial PointSource™ Optics: This proprietary technology makes it easy to enter undilated pupils, and provides a 25° field of view resulting in a view of the fundus that's 5X larger than you see with a standard ophthalmoscope in an undilated eye.

Patented Glare Extinguishment System: Prevents interference from unwanted glare and reflection.

Dynamic Focusing Wheel: Adjusts the focus in a continuous, smooth action for more precise control and optimum view. Available focusing range is -20 to +20 diopters.

3.5 V Lithium Ion Power Handle:

The 3.5 V Lithium Ion Power Handle features “Smarter” Lithium Ion battery technology that powers any 3.5 V instrument head for 110 minutes—twice as long as the conventional handle. Its rugged plastic resin construction and flexible, ergonomic design make the 3.5 V Lithium Ion Power Handle more comfortable for any sized hand.

- Half the weight of current 3.5 V handles
- Average 110 minutes of on-time between charging
- Smallest diameter of any full-size 3.5 V power handle
- Low-battery indicator
- No need to discharge between charging cycles

The Panoptic ophthalmoscope also includes a charger and zip case.

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	WA11824-VSM-cl
Diagnostic Head Type	Ophthalmoscope

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights