

## REUSABLE WARMING BLANKETS

CWS-VET0-CONFIG\*



### PRODUCT DESCRIPTION

---

## Reuseable warming blanket made from easy-to-clean, high-quality polyester fabric

Reuseable warming blanket made from soft touch, high-quality polyester. A 5" diameter opening with a nylon cord secures the blanket to all major brands of warm air convection units. The forced air warmers significantly reduce the risk of patient hypothermia during and after anaesthesia procedures.

The inlay is a durable and waterproof nylon material which circulates the airflow upwards towards the patient laying down. This durability comes from the thickness of the material and the strength of the stitching. The design of the blanket makes it safe and easy to clean. Unlike most warming blankets, these are machine washable (hang to dry) and more environmentally conscious than their disposable counterparts.

**Small 30"L x 20"W**

**Medium 36"L x 24"W**

**Large 48"L x 24"W**

**Small 75cm L x 70 cm W**

**Medium 90cm L x 60 cm W**

**Large 120cm L x 60 cm W**

Choose from Small, Medium or Large from the range of reusable warming blankets.

#### Key Features:

- 3 Available sizes for different shaped and size mammals
- Thick, durable, high-quality waterproof warming blankets
- Retains heat and warm to the touch
- Soft, comfortable fabric ideal for using under patients pre or post-op.
- Heavy-duty nylon and polyester is easy to clean and resists damage/punctures
- Reusable - Economical and saves money
- Hang dry (Not suitable for tumble-dryers)
- Even distribution of heat
- Fast warm up time - forced air warmers that get to temperature in minutes
- Designed to minimise excess air being blown around the surgical site.

## ADDITIONAL INFORMATION

<b>SKU</b>	CWS-VETO-config*
------------	------------------

<b>Animal</b>	Small Animal
---------------	--------------

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights