

MINDRAY SPO₂ VET SENSOR CLIPS

IMEC-115-0-GROUPED



PRODUCT DESCRIPTION

Mindray SpO₂ Vet Sensor Clips

Extra or replacement clips for your SpO₂ sensors, helping you to avoid replacing the entire sensor in the event of damage. Additional clips can also facilitate easier changes between patients. Compatible with all Mindray SpO₂ sensors, these clips are suitable for connecting to a patient's tongue, ears, or between pads.

Two sizes available:

- Small = 4.5 cm
- Large = 6 cm

ADDITIONAL INFORMATION

SKU	IMEC-115-0-grouped
Accessory Parameter	SpO ₂

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights