

ADAPTABLE ANAESTHESIA MACHINE MECHANICAL VENTILATION UPGRADE

102-068-WEBINAR-PROMO



PRODUCT DESCRIPTION

An update on a classic from the Burton's Adaptable range, the Merlin Shelf helps your practice integrate ventilation.

Built upon the original Burton's AAA (Adaptable Anaesthesia Apparatus) frame concept, the Merlin shelf is designed to go on any Burton's AAA machine to improve the adaptability and versatility of the frame.

The Anaesthesia machine sits on top, leaving enough room to fit the shelf underneath to house a Merlin Small Animal Ventilator, with the potential for an Oxygen Concentrator to go underneath. The 16Kg Merlin can be supported by this robust purpose-built shelving unit, which is made from easy to clean, high-grade 304 stainless-steel. The shelf can also be used to house and store other essential items.

This shelf design has a small footprint, which fits perfectly into the existing compact footprint of the AAA machine frame, enabling your practice to put together the perfect bespoke unit configuration.

Bundle Includes:

- **Intergrated Merlin Shelf**
- **Vetronic Merlin Ventilator**
- **Ventilation Webinar for your practice**

Key Features:

- Complete Burton's Lifetime stainless-steel frame – durable, robust, easy-clean.
- Designed and manufactured by Burton's in the UK.
- Standard 1 year guarantee.

- Anaesthesia Machine, Concentrator and vaporiser(s) sold separately.

Technical Specification:

Dimensions: Height 251mm X Width 455mm X Depth 316mm

Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent. Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

Vaporiser & AAA Sold Separately

ADDITIONAL INFORMATION

Animal	Mixed
SKU	102-068-webinar-promo
Catagory Sub Heading	Includes Ventilator shelf, Vetronic Merlin Ventilator & Clinical ventilation webinar
Brand	Burtons
Gas Type	Various Options
Gas Supply	'F' Size Cylinder(s)
Flowmeter Type	Flowmeter
Concentrator	No
Optional O2 Alarm	No
Optional N2O	No
Emergency O2 Flush	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

