

## SPITZ FULL HD CAMERA SET

END-300-001-480



### PRODUCT DESCRIPTION

---

## Introducing the revolutionary SPiTZ range, offering comprehensive solutions for veterinary diagnostics and procedures.

Mobile endoscopy becomes possible with this High Definition camera that allows endoscopic examination with all required parts in one complete set. Integrated illumination and a high-performance camera head, powered by an external Desktop-Computer or Laptop via USB 3.0 makes examination and clinical surgery possible.

#### Set includes:

- Camera head
- Endoscope coupler
- Light Source
- Suitcase
- Software
- TV-adapter

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	END-300-001-480

<b>Brand</b>	SPITZ
<b>Procedure Type</b>	Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Endoscopy Device</b>	USB Camera system
<b>Animal Type</b>	Companion animals, Equine and Large Animal
<b>Image Resolution</b>	Full HD

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights