

SPITZ LUMEN LED 1 LIGHT SOURCE SYSTEM

END-300-900-180



PRODUCT DESCRIPTION

SPiTZ Lumen LED1 Light Source System

Key Features:

- Daylight intensity display
- 20.000h lamp lifetime
- Manual or automatic light intensity control
- Waterproof touch control panel
- Low heat conduction by improved LED technology

System includes:

- Main Unit
- BNC cable
- Mains cable
- Operating Manual

ADDITIONAL INFORMATION

Animal	Companion animals, Equine and Large Animal
SKU	END-300-900-180
Brand	SPITZ

Procedure Type	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
Endoscopy Device	Light source
Animal Type	Companion animals, Equine and Large Animal
Image Resolution	LED

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights