

## SPITZ ADAPTER FOR F/O LIGHT CARRIER FOR TYPE LUXTEC LIGHT SOURCE

END-250-048-000



### PRODUCT DESCRIPTION

SPITZ Adapter for F/O Light carrier for Type Luxtec Light source

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	END-250-048-000
<b>Brand</b>	SPITZ
<b>Procedure Type</b>	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Endoscopy Device</b>	Adaptor piece, light guide cable
<b>Animal Type</b>	Companion animals, Equine and Large Animal

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

