

SPITZ ARTHRO P 1 IRRIGATION PUMP

END-225-906-000



PRODUCT DESCRIPTION

SPITZ ARTHRO P 1 Irrigation Pump For Arthroscopy

Key features:

- Hydrostatic pressure regulation Silicone Tubing Set included (END-225-910-001)
- Easy to operate, splash-proofed display
- Highly reliable tube guide and handling
- Skid-proof pedestals

Unit contains:

- Main Unit
- Tubing Set
- Mains Cable
- 2 Spare Fuses
- Operating Manual

ADDITIONAL INFORMATION

Animal	Companion animals, Equine and Large Animal
SKU	END-225-906-000
Brand	SPITZ
Procedure Type	Arthroscopy
Endoscopy Device	Arthroscopy device
Animal Type	Companion animals, Equine and Large Animal

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights