

SPITZ STERILE FILTER WITH LL-CONNECTOR

END-300-902-007



PRODUCT DESCRIPTION

SPITZ Sterile Filter with LL-Connector package of 50pcs.

ADDITIONAL INFORMATION

Animal	Companion animals, Equine and Large Animal
SKU	END-300-902-007
Brand	SPITZ
Procedure Type	Minimally invasive surgery, Laparoscopy
Endoscopy Device	Laparoscopy Device, Insufflator
Animal Type	Companion animals, Equine and Large Animal

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

