

SPITZ VETERINARY OTOSCOPE 5 CHARR, 4MM, 65MM WORKING LENGTH

END-410-005-000



PRODUCT DESCRIPTION

SPiTZ fully autoclavable Veterinary Otoscope

- Working Length: 65mm
- Diameter: 4mm
- Angle of View: 0°
- 5 Charr

ADDITIONAL INFORMATION

Animal	Companion Animals, Exotics
SKU	END-410-005-000
Brand	SPiTZ

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

