

SPITZ SINUSCOPE/ RHINOSCOPE 30°, 2.7MM, 180MM WORKING LENGTH

END-252-127-030



PRODUCT DESCRIPTION

SPITZ fully autoclavable Veterinary Sinuscope/ Rhinoscope - Working Length 180mm - Diameter: 2.7mm - Angle of view: 30°

ADDITIONAL INFORMATION

| | |
|---------------|----------------------------|
| Animal | Companion Animals, Exotics |
| SKU | END-252-127-030 |
| Brand | SPITZ |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights