

## SPITZ RZ CYSTOSCOPE, TRANSURETHRAL INJECTION, 12° 4.0MM DIAMETER, 300MM WORKING LENGTH

END-253-804-012



### PRODUCT DESCRIPTION

## SPITZ fully autoclavable Veterinary Cystoscope for transurethral Injection

- Working Length : 300mm
- Diameter: 4.0mm
- Angle of view: 12°

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals
<b>SKU</b>	END-253-804-012
<b>Brand</b>	SPITZ

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

