

**SPITZ ADVANCED SEMI-RIGID URETERO-RENSCOPE, 425MM  
WORKING LENGTH**

END-253-915-540

**PRODUCT DESCRIPTION****SPiTZ fully autoclavable Spitz Advanced Semi-Rigid Uretero-  
Renscope for artificial insemination**

- Disital Tip: 4.5 Charr.
- Distal: 5.0 Charr.
- Proximal: 6.5 Charr.
- Working Chanel 4.0 Charr.
- Field of View 65°
- Direction of View 5°

**ADDITIONAL INFORMATION**

<b>Animal</b>	Companion animals
<b>SKU</b>	END-253-915-540
<b>Brand</b>	SPiTZ

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights