

## LAPAROSCOPE 4K 0°

END-300-155-400



### PRODUCT DESCRIPTION

SPITZ fully autoclavable 4k Laparoscope - Working Length : 300mm - Diameter: 5mm - Angle of view: 0°

### ADDITIONAL INFORMATION

|               |                                   |
|---------------|-----------------------------------|
| <b>Animal</b> | Companion Animals,Equine, Exotics |
| <b>SKU</b>    | END-300-155-400                   |
| <b>Brand</b>  | SPITZ                             |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights