

WELCH ALLYN OPERATING LITHIUM ION RECHARGEABLE DIAGNOSTIC SET

WA96274-VSM



PRODUCT DESCRIPTION

The 719 Lithium Power handle features smart lithium ion battery technology, which is environmentally friendly and has an average usage time of 120 minutes

Coaxial Ophthalmoscope:

- 6 apertures for all general and specialist use
- 28 lenses -25 to +40 diopters for better resolution
- Rubber brow rest prevents scratching of eyeglasses
- Illuminated lens dial clearly identifies settings
- Coaxial optical system for shadow-free spot and easier entry into undilated pupil
- Crossed linear polarizing filter to reduce corneal reflection
- Red-free filter for use with any aperture
- 18 possible aperture combinations for greater versatility
- Sealed optics for better dust protection

Operating Otoscope:

- Open system for convenience during procedures
- Rotatable lens and speculum for ease of use

Fiber optic illumination for cool light

3.5 V Lithium Ion Power Handle:

The 3.5 V Lithium Ion Power Handle features “Smarter” Lithium Ion battery technology that powers any 3.5 V instrument head for 110 minutes twice as long as the conventional handle. Its rugged plastic resin construction and flexible, ergonomic design make the 3.5 V Lithium Ion Power Handle more comfortable for any sized hand.

Half the weight of current 3.5 V handles

Average 110 minutes of on-time between charging*

Smallest diameter of any full-size 3.5 V power handle

Low-battery indicator

No need to discharge between charging cycles

The operating lithium-ion rechargeable diagnostic Set also includes 3 x polypropylene specula and (4,5 & 7mm) as well as a zip case, with the charger and charging pod included in the set.

ADDITIONAL INFORMATION

Animal	N/A
SKU	WA96274-VSM
Generic Size	X-Small
Digital	No
Diagnostic Head Type	Operating Otoscope
Battery Power Type	Lithium Rechargeable

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights