

BR-M25 BRONCHOSCOPE

END-BR-M25



PRODUCT DESCRIPTION

The BR-M25 Semi-Disposable Flexible Bronchoscope is a state-of-the-art tool engineered for the diagnosis and treatment of respiratory conditions in veterinary patients.

Semi-disposable scopes are a cost-effective solution to limiting contamination risks, incurring less processing costs, and reducing time spent on maintenance whilst remaining the ideal solution for procedures. Each scope has up to 10 uses.

Combining advanced design with high-quality imaging, it is perfectly suited to veterinary use, delivering exceptional precision in airway examinations, biopsies, foreign body retrieval treatments, and surgeries. This versatile bronchoscope is specifically designed for procedures involving the nasal cavity, upper respiratory tract, trachea, and bronchus, ensuring reliable performance across a range of diagnostic and therapeutic applications.

It has a working length of 600mm with 2-way rotation making it the ideal choice for veterinarians. The insertion tube outer diameter for the endoscope is 2.5mm, with a minimum working channel inner diameter of 1.2mm. The imaging quality is at HD resolution (400 x 400 pixels) for crystal-clear visuals to help you navigate challenging procedures. The field of view for the scope is $120^{\circ} \pm 10\%$, with a $0^{\circ} \pm 10^{\circ}$ forward-facing direction of view and a depth of field ranging from 3mm to 50mm.

The BR-M25 provides high-resolution imaging via the MS-8 Video Processor, ensuring clear, detailed visuals for examinations of the nasal cavity, upper respiratory tract, trachea, and bronchus. The unit can be used on small animal, equine, wildlife, zoo, avian and exotics. It can be utilised in a mixed practice or even in a remote setting.

The lightweight handle of the unit includes a suction port for seamless integration with irrigation systems and a biopsy port for channel instruments, with 2 customisable buttons that can be programmed from the all-in-one MS-8 processor to control photographing, videoing, white balance adjustments and graphics freezing with the selected buttons displayed as hotkeys.

There is a manual lever at the base of the handle which controls the bending section and the distal tip of the scope for precise

navigation. The high-quality medical grade patient tube is equipped with an integrated camera and LED light source. These components deliver high-definition visuals, transmitting image data to the MS-8 Image Processor via a 2500mm connecting wire.

Every bronchoscope is packed in sterilised, single-use packaging with the date of manufacture and expiration date clearly indicated on the label, ensuring an aseptic shelf life of 3 years. This minimises cross-contamination risks, enhancing safety for both patient and practitioner. Professional disinfectant wipes that are effective against bacteria, yeasts, viruses, and fungi which are tested to veterinary standards (European PT3 test standards for microbiological testing on key veterinary equipment) are suitable for cleaning the outside surfaces of these scopes.

With seven models within the BR series, there is ample choice to meet the unique and diverse needs of veterinary practitioners who are making diagnosis and undertaking treatments on patients. The various sizes satisfy the requirements of different operational settings. This range is best suited to alveolar lavage, airway examinations, biopsies, foreign body removal, intubation, medication delivery and secretion suction procedures. The units are fast to set up, with the ability to run on battery for over 4 hours, making them the ideal choice for out of hours and remote work.

The BR-M25 bronchoscope offers reliability, flexibility, and high performance, making it a trusted choice for veterinary professionals seeking excellence in endoscopic care. It is compatible with all the conventional accessories such as biopsy forceps, foreign body forceps, baskets, prong graspers, brushes, and mucus collection traps. and 1.0mm instruments used in endoscopic care.

Key Features

- HD Resolution for clear, high-quality imaging.
- Built-in LED light source and high-definition camera.
- Variable field of view and sharp focus in narrow spaces and at great depths.
- High-quality patient tube with flexible distal end for smooth navigation.
- Supports a wide range of procedures versatile applications and advanced techniques.
- Easy and wide aspiration and sampling through the working channel.
- User-friendly multi-functional handle with programmable hotkeys for ease of use.
- Manual lever for precise control and adjustment when using the bending section.
- Safe and hygienic sterilised packaging with a 3-year aseptic shelf life, reducing cross-contamination risks.

Specification

Endoscope Model	BR-M25 - Semi-Disposable Flexible Video Bronchoscope
Field of View (°)	120° ±10%
Direction of View (°)	0°±10°
Depth of Field (mm)	3mm - 50mm
Insertion tube diameter(mm)	2.5mm
Working Length (mm)	600mm±10%
Working Channel (mm)	1.2mm±10%
Bendable Portion (°)	210° (upwards) and 210° (downwards), with a deviation of -5% (excluding upper limit)
Functions	Photographing / Videoing / White Balance / Graphics Freezing
Connections	USB / HDMI
U/D, L/R rotation(ways)	2 Way Rotation
CMOS Pixel (px)	400px*400px
Optical	HD, Wide Angle
Overall Length (mm)	860mm±10%
Operational Environment Temperature (°C)	5-40°C(41-104°F)

Supporting Products:

- END-310-901-120 - Biopsy forceps Long Oval Cup, fenestrated D: 1.0 mm, 120 cm length
- END-310-960-120 - Alligator Forceps flex. foreign body, Forceps D: 1 mm, WL 120 cm
- END-NB3FR120 - 1.0mm x 1200mm Basket
- END-SG3FR120 - 1.0mm x 1200mm 3 Prong Grasper
- END-RC3FR120-S - 1.0mm x 1200mm pack of 10 sterile BAL catheters
- END-253-962-003 - Brush for Cytology, 3 Charr., For single use
- END-310-901-120 - Biopsy forceps Long Oval Cup, fenestrated D: 1.0 mm, 120 cm length
- END-17021182 - Mucus Collection Trap 30ml - Pk of 10
- END-207230420 - Cleaning Brush Ø2mm, length 2300mm, Disposable Double-Sided, Ball Head, Tube Ø1.0mm

List of Other Semi-Disposable BR Scopes

Semi-disposable BR scopes							
Model	BR-M58H	BR-M58H-1500	BR-M52H	BR-M52H-1000	BR-M40H	BR-M32	BR-M25
Insertion tube diameter(mm)	6.0	6.0	5.0	5.0	4.0	3.1	2.5
Working channel(mm)	2.8	2.8	2.2	2.2	1.2	1.2	1.2
Working length(mm)	600	1500	600	1000	600	600	600
U/D, L/R rotation(ways)	2	4	2	4	2	2	2
CMOS Pixel	1280*720 HD	1280*720 HD	1280*720 HD	1280*720 HD	1280*720 HD	400*400	400*400

ADDITIONAL INFORMATION

Animal	Companion animals
SKU	END-BR-M25
Overall Product Dimensions (mm)	860
Colour/Style	White
Procedure Type	Bronchoscopy, Diagnostic Endoscopy , Respiratory Examination, Rhinoscopy, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
Biopsy Channel	1.2mm
Working Length (mm)	600
Fitting	Universal

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights