

**BURTONS M880B HANDHELD CAPNOGRAPH AND PULSE OXIMETER
- CLEARANCE**

129-870-CL-2-2

**PRODUCT DESCRIPTION****Ultra clear, ultra lightweight touchscreen Capnograph & Pulse Oximetry monitor. The perfect addition to any practice!**

Boasting a large 4.3inch, high resolution, auto- rotation touchscreen and a large capacity lithium battery. The M880B is designed for continuous monitoring or spot-checking of CO2 and SPO2 signals in patients.

Key Features:

- Extremely light, for carrying and handling Only 0.5KG with full configuration, including battery.
- 4.3inch TFT high resolution with touch screen & auto-rotation for user ease of use with any orientation.
- Protection against ingress of liquid (IPX2 Waterproof)
- Specially designed water filter, ensuring humid environments don't effect readings
- Configurable for flux for horse, dog and cat
- Large capacity lithium battery for up to 18 hours use
- CO2 start/off button to save power
- Selectable layout modes including font and waveforms
- Visual and audio alarm for exceeding limits and Apnea; pitch tone

Size:	4.3inch
Resolution:	480*272
Touch:	Resistive Touch
Auto-rotation Direction	4 Direction

Power Adaptor:	Input: AC 100 V±2 240V Output: DC 5V/2A
CO2:	
Measurement Way:	Infrared Spectrum
Measurement Range:	0-19.7% (0-150mmHg 0-20kPA)
CO2 Accuracy:	<5.0% CO2 ±2mmHg >5.0% CO2 <6% of the reading
CO2 Resolution:	0.1mmHg
Gas Flow Rate:	60-80ml/min
Response Time:	1/4±3s
Calibration:	Offset calibration, auto/ manual; gain calibration
SpO2	
Measurement Range:	0-100%
Resolution:	1%
Accuracy:	70-100%± 2% 35-69%± 3% 0-34% definition

S/N: M67B0001104

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	129-870-cl-2-2
Colour/Style	White
Brand	BLT
Monitor Parameters	CO ₂ , SpO ₂
Screen Type	Touchscreen
Power Source	Battery
Optional Modules	N/A

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights