

ASTRAMAX HD MOBILE MOUNTED OPERATING LIGHT - CLEARANCE

AMAX30FM-CL



PRODUCT DESCRIPTION

The aim of the development of Astramax was to bridge the gap between the biggest operating theatre lights and smaller minor operating lights.

Astramax carries all of the advantages of HD-LED technology including:

Full Spectrum Colour Rendition

The first LED lighting with near perfect colour rendition across the full visible spectrum. This reduces eye strain and enables small differences in tissue to be distinguished easily.

Fat Beam Illumination

More light flux provides the biggest beams of light in the business, bigger illuminated fields and more light across the full width of the illuminated area.

Red Balance Control

The highest R9 red colour rendition of any product provides the optimum visualisation of red tissues. Red balance enhancement compensates for our natural weakness in distinguishing shades of red.

- HD-LED Technology
- High performance compact surgical light
- LED life of 40,000 hours
- 140,000 Lux
- Large movement range
- Economical
- Environmentally friendly
- Lamp head sealed to IP54
- Shadow control
- 5-100% light intensity control

Data Common To All Models

| | | |
|--------------------------------|------------------|---------------|
| Lux | 1m | 140,000 |
| Colour Rendition Index | Ra | 95 |
| Colour Rendition (Red) | R9 | 98 |
| Colour Temperature | Kelvin | 3,700 - 5,000 |
| LED Life | Hours | ≥40,000 Hours |
| Power Consumption | Watts | <50 |
| Endoscope Settings | | 5% |
| Beam Size | D10 | 140-310mm |
| Weight (Ceiling, Wall, Mobile) | 15Kg, 11Kg, 30Kg | |

ADDITIONAL INFORMATION

| | |
|------------------|-------------|
| Animal | N/A |
| SKU | AMAX30FM-cl |
| Head Type | Single |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights