

**MACH LED CEILING LIGHT - FOCUSABLE - STERILISABLE HANDLE - 25CM - CLEARANCE**

MACHLED130F-C-SH-25-CL

**PRODUCT DESCRIPTION**

MACH LED Ceiling Light - Focusable - Sterilisable Handle - 25cm

MACH LED 130F Ceiling Mounted Light  
with 25cm Ceiling Pole Kit  
Includes Removeable Sterilizable Handle

\*\*\*Please note: Anchorage Ring is not included\*\*\*

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	MACHLED130F-C-SH-25-cl
<b>Brand</b>	Dr Mach
<b>Autoclavable Handle</b>	Yes
<b>Focusable</b>	Yes

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights