

HUGER 40W LIGHT SOURCE WITH HARD CARRY CASE - CLEARANCE

END-HGVIS-68-MON-PX-CL



PRODUCT DESCRIPTION

The VIS-68 endoscopy processor in combination with HQ scopes delivers the sophisticated endoscopic images via improved image processing technology, enhanced colour reproduction, minimized image noise and reduced halation

Part Exchange HUGER 40 W LED Light Source
with 12" Intergrated IPS Monitor
No Warranty

Huger 40 W LED Light Source with 12" Intergrated IPS Monitor. Hard carry case (END-FC03) included.

IPS display screen

The IPS display screen with 60fp/s capacity presents the optimal colour reproduction for both still and dynamic endoscopic images

Power Supply:

- Voltage: 100-240V AC
- Frequency: 50~60Hz
- Consumption electric power: 150VA

Size:

- Dimensions (W x H x D): 375*320*130mm
- Weight: 9.5kg(with monitor)

Observation:

- Digital signal output: DVI*1, HDMI*1
- White balance adjustment: Automatic WB
- Colour adjustment: ± 5 , total 11 levels
- Colour saturation adjustment: 0~10, total 11 levels
- Brightness adjustment: 0~10, total 11 levels

Documentation:

- Image capture: Yes
- Video record: Yes
- Image and video storage: SD card compatible
- Data management: Assign one file for each case

Illumination:

- Examination lamp: LED lamp 40W
- Average lamp life:35000 hours
- Colour temperature: 5300K

S/N: VHD681017

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-HGVIS-68-MON-PX-cl
Brand	Huger
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Illumination Type	L.E.D

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights