

## MACH LED 130 CEILING MOUNTED LIGHT WITH 60CM CEILING KIT WITH STERISABLE HANDLE- CLEARANCE

MACHLED130-C/KIT-60-SH-CL



### PRODUCT DESCRIPTION

---

## MACH LED 130 Focusable Sterisable Handle Ceiling Light - With Lamp Head - 60cm

Mach LED 130 Focusable Light & Ceiling Kit - Pole length set at 60cm

Consists of:

- Ceiling Tube
- Ceiling Flange
- Swing Arm
- Sterilisable Handle
- Mach LED 130 Lamp Head For Ceiling Mounted Light

**Approx ceiling height of 2.8m**

\*Please note: Anchorage Ring is not included

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	MACHLED130-C/KIT-60-SH-cl
<b>Brand</b>	Dr Mach

**Autoclavable Handle**

Yes

**Focusable**

No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights