


## BURTONS SUPERBLUE GIVING SETS - STANDARD NEEDLE FREE VALVE

129-012\*BUNDLE



- 1 Drip chamber
- 2 Translucent blue tubing for pump placement
- 3 'C' Clip flow shut-off clamp
- 4 Needle-free valve
- 5 Translucent tubing
- 6 Blue precision 'V' roller clamp
- 7 Translucent tubing
- 8 360 spin rotating luer lock
- 9 Prime Stop Cap

### PRODUCT DESCRIPTION

## Burtons SuperBlue Giving Sets - Standard Needle Free Valve

Our high-quality Burton's SuperBlue Giving Set as part of our Infusion Perfection range - A standard needle free valve design exclusive to Burton's, the needle free valves are compatible with both standard syringes and syringes with a luer lock for a secure fit. This SuperBlue giving set has been rigorously tested by our experts to ensure it meets the highest standards of quality and performance. Its innovative design includes a blue placement indicator colour coding. The ease of use makes these giving sets the standout choice for veterinary professionals. They are ideal for short term infusions (1-2 days) and can be used with all infusion pumps as well as via a gravity feed. Each giving set is single use only.

### Key Features:

- **Blue Placement Indicator:** A unique, blue-tinted section on the line ensures easy identification of the correct placement for infusion pumps, streamlining setup and reducing errors. Blue pump section is 900mm.
- **Compatibility:** Specifically designed to work seamlessly with the Burton's range of infusion pumps, including the EN-V5, EN-V5+, IVY, Niki IV, as well as other infusion pumps. The set can also work when gravity-fed.
- **Ultrasonic Sensor-Friendly:** Custom-built to avoid interference with ultrasonic sensors on pumps, ensuring smooth operation and preventing occlusions.
- **Instant Flow Control:** Equipped with a 'C' Clip flow shut-off clamp for immediate, temporary flow cessation when needed.
- **Precision 'V' Roller Clamp:** Allows for accurate flow rate adjustment.
- **Anti-Kink Resistant Tubing:** Ensures consistent fluid delivery and seeks to further prevent occlusions, even with active patients.
- **Durable High-Quality Materials:** Made from DEHP-free PVC tubing, latex-free components, and high-quality materials such

as PE, PTFE, ABS, and nylon for reliable performance.

- **Rotating Luer Lock:** Enhances versatility and ease of use during procedures. The rotating luer lock also allows for easier connection to IV cannulas.
- **Needle-Free Valves:** Includes two clear blue needle-free valves for added safety and convenience.

## Why Choose Burtons SuperBlue Giving Sets?

Exclusive to Burtons, this giving set has been rigorously tested by our experts to ensure it meets the highest standards of quality and performance. Its innovative design, including the blue placement indicator and anti-kink tubing, makes it a standout choice for veterinary professionals. All Burtons sets are designed in such a way to ensure consistent calibration across the range, allowing a single calibration value – There is no need to remember to change your settings between using Burtons sets!

## Specifications:

- **Drip Chamber:** 20 drops/ml, 15µm (Micron) filter
- **Clamps:** Blue POM clamp, blue precision 'V' roller clamp (ABS), 'C' Clip flow shut-off clamp
- **Needle-Free Valves:** 2 clear blue PC/silicone valves (latex-free)
- **Prime Stop Cap:** Blue PP/PTFE
- **Priming Volume:** 13ml
- **Patient Connector:** Transparent male luer lock (ABS)
- **Overall Length:** 2.07m
- **Dropper:** Approx 20 drops/ml
- **Solution Filter:** 15µm (Micron)
- **Sterilised:** Sterilised with ETO and given a dated two-year shelf life from the manufacture date.

## Enhance Your Equipment:

Ensure accurate readings by having your infusion pump calibrated by our expert workshop team. Contact us today to discuss your needs!

Burtons Giving Sets Comparison Chart				
Part number	129-010	129-011	129-012	129-013
Name	Coiled	Paediatric	Needle free	Standard
Total length	4.99m	1.71m	2.07m	2.91m
Pump section	✓	✓	✓	✓
Injection ports	2	1	✗	2
Needle-free valves	✗	✗	2	✗
Prime volume	24.8ml	10.8ml	13ml	19.4ml
C clamp	✓	✓	✓	✓
V roller	✓	✓	✓	✓
Drops/ml	20	60	20	20
Uses	Large patients and walk-in kennels	Small patients	Prevent needle stick injuries	Everyday

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	129-012*bundle
<b>Brand</b>	Burtons
<b>Category Sub Heading</b>	Save 20% on packs of 100!

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights