

EXAMION FIRECR DENTAL SYSTEM - WITH FIRECR READER, GENERATOR + AQS SOFTWARE

XR-EXYFIRE3-PROMO25-BUNDLE



PRODUCT DESCRIPTION

This fully featured CR dental system delivers high-quality dental images quickly and affordably for the veterinary clinic.

This bundle includes The Examion FireCR Reader + Examion AQS Software that can be added to your existing PC + a wall-mounted generator with a standard arm.

Upgrade to the Premium Fire CR Reader for larger plate support, enabling rapid examination and reduced radiation exposure. The large size 4 plate can also be used to image extremities.

This dental CR solution has been designed to use with your existing practice PC and includes:

- FireCR reader
- Imaging plate kits of sizes 4 x 2 & 1 x 4c with 100 protective sleeve sizes for both sizes.
- Cleaning fluid
- Examion AQS software including the following modules: X-ray, FireCR control, Image optimisation, 1 x license, AQS PACS.
- Static Dental Wall Mounted Generator with standard arm (184cm)

The FireCR's compact size makes it suitable in even the smallest examination rooms, without compromising quality. It rapidly reads reusable, easy-to-position intraoral imaging plates, generating high-quality X-ray image results.

Working alongside the included dental x-ray generator, the reader's built-in FireCR software includes bespoke processing and archiving modules, allowing you to simply purchase the Fire CR reader to complete your dental setup.

Imaging plates are processed by being placed directly onto the feeder, which magnetically holds them in place. This contact-free reading ensures a longer plate life and no additional cassette adapters are needed to read the imaging plate.

Key Features:

- Simple, two-step process – after exposure, the Fire CR provides fast, reliable images.
- Dental image plates are available in different dental sizes for different animal species.

- It's rapid, high-quality reading of re-usable, easy-to-position bitewing, intraoral positioning plates instantly enhances your dental x-raying.
- Exceptionally compact at just 12cm wide and 26cm high, suitable for the smallest of spaces

Specifications:

- Supported Intraoral Imaging Plate sizes: Sizes 0 - 4c
- Pixel size: 35µm
- Sampling Pixel Pitch: 14.3 lp/mm
- Readout Time: 4.1 - 7.2s
- Resolution: Data Capture: 16bits/pixel, 65,000 graytones
- Eraser: Embedded
- Dimensions: 120mmW x 265mmH x 318mmD
- Weight: 4.7kg
- System configuration: table top
- Computer workstation minimum requirements: Intel Dual Core processor / 2GB RAM / 10GB free disc space/ Windows 7,8,10

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	XR-EXYFIRE3-PROMO25-BUNDLE
Generic Size	Small
Brand	Examion
Catagory Sub Heading	Dental CR System ready to add to your PC

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights