

OPTICLAR DESK DIAGNOSTIC SET WITH OPHTHALMOSCOPE & OPERATING OTOSCOPE

OPT100.020.034



PRODUCT DESCRIPTION

Veterinary Desk Ophthalmoscope & Operating Otoscope Diagnostic Set

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

Set includes:

- Operating Otoscope with 3.0x magnification
- L28 Practice ophthalmoscope
- 2 Adapt rechargeable handle
- Twin port charger
- 15 plastic Specula tips (5 of each size) - 57.0mm x 4.0mm diameter, 65mm x 6.0mm diameter, 90.0mm x 7.0mm diameter
- Zip case

Key Features

- TrueTone LED
- Excellent LED illumination

- Bayonet tip connection
- OpenHeadDesign - allows introduction of instruments without obscuring the view and for easier instrument access.
- Distal fibre optics gives cool light to examination site with no reflections and provides excellent view of the tympanic membrane
- Swivel viewing lens cannot be lost and is guarded on lower edge from potential breakage
- Wide angle viewing window for easy instrumentation under magnification
- Head available with German or American tip fitting
- 28 lens from -20d to +40d
- Durable, solid metal construction
- Opticlar and Welch Allyn connection
- Compatibility with Heine veterinary tips (Heine Beta 200 Veterinary style)
- 5 apertures: Large, medium and small spot; Graticule and Slit
- 10 year LED warranty
- 20,000 hour plus LED life guarantees unlimited use with no blown bulb downtime
- 5 year product warranty
- Supplied with 3 plastic tips

Adapt Desk Rechargeable Handle

- Lithium Ion
- Robust metal sleeve
- Knurled surface for maximum grip
- Rotary control allows the LED brightness to be easily adjusted

Recharging Pod

- Comply with national and international standards
- CE marked
- Single port recharging pod
- Supplied as standard with 3 pin UK plug, other configurations available upon request
- Generous 3m/ 10ft lead

Specifications:

- Charging time: 12 hours via usb
- Running duration: > 120 hours

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT100.020.034
Brand	Opticlar
Catagory Sub Heading	Save 16%
Diagnostic Head Type	Co-axial Ophthalmoscope & Operating Otoscope

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights