

OPTICLAR DIAGNOSTIC CASE SET - WITH OPHTHALMOSCOPE AND SLIT OTOSCOPE

OPT700.020.020D



PRODUCT DESCRIPTION

Veterinary Diagnostic Set - C-Cell Battery Handle

Introducing a new standard in LED Veterinary Diagnostic Equipment – The most popular and widely used veterinary diagnostic set in the Opticlar range.

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

Set includes:

- S1 Veterinary LED slit otoscope
- L28 Practice ophthalmoscope with cobalt blue filter, light sensitive filter and 5 apertures
- 1 Adapt C-Cell handle
- 15 plastic Specula tips (5 of each Size) - 57.0mm x 4.0mm diameter, 65mm x 6.0mm diameter, 90.0mm x 7.0mm diameter
- Zip case

Key Features

- TrueTone multi-voltage LED slit head
- Can be 'locked' on to handle

- Excellent LED illumination
- Bayonet tip connection
- Swivelling 3x magnification lens
- Durable, solid metal construction
- Opticlar and Welch Allyn connection
- 10 year LED warranty
- 5 year product warranty
- Supplied 15 plastic tips (5 of each Size) - 57.0mm x 4.0mm diameter, 65mm x 6.0mm diameter, 90.0mm x 7.0mm diameter

Specifications:

- Running duration: 125 hours

C-Cell Battery Handle

- C-Cell Battery Powered
- Robust metal sleeve
- Knurled surface for maximum grip
- Rotary control allows the LED brightness to be easily adjusted

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.020.020D
Brand	Opticlar
Category Sub Heading	Save 15%
Diagnostic Head Type	Co-axial Ophthalmoscope, Slit Ootoscope

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights