

## OPTICLAR DIAGNOSTIC DESK SET - WITH OPHTHALMOSCOPE AND SLIT OTOSCOPE - SINGLE DOCK

OPT700.020.022



### PRODUCT DESCRIPTION

---

## Veterinary Diagnostic Set - Single Port Desk Charger

Introducing a new standard in LED Veterinary Diagnostic Equipment - The most popular and widely used veterinary diagnostic set in the Opticlar range.

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

#### Set includes:

- S1 Veterinary LED slit otoscope
- L28 Practice ophthalmoscope with cobalt blue filter, light sensitive filter and 5 apertures
- 1 Adapt rechargeable handle
- Single port charger
- Supplied with 3 metal tips (1 of each size) - 57.0mm x 4.0mm diameter, 65mm x 6.0mm diameter, 90.0mm x 7.0mm diameter
- Zip case

#### Key Features

- TrueTone multi-voltage LED slit head

- Can be 'locked' on to handle
- Excellent LED illumination
- Bayonet tip connection
- Swivelling 3x magnification lens
- Durable, solid metal construction
- Opticlar and Welch Allyn connection
- 10 year LED warranty
- 5 year product warranty
- Supplied with 3 metal tips

### Adapt Desk Rechargeable Handle

- Lithium Ion
- Robust metal sleeve
- Knurled surface for maximum grip
- Rotary control allows the LED brightness to be easily adjusted

### Recharging Pod

- Comply with national and international standards
- CE marked
- Single port recharging pod
- Supplied as standard with 3 pin UK plug, other configurations available upon request
- Generous 3m/ 10ft lead

### Specifications:

- Charging time: 2 hours
- Running duration: > 20 hours

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	OPT700.020.022
<b>Brand</b>	Opticlar
<b>Catagory Sub Heading</b>	Save 16%
<b>Diagnostic Head Type</b>	Co-axial Ophthalmoscope, Slit Oscope

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

