

OPTICLAR DIAGNOSTIC BATTERY CASE SET WITH OPHTHALMOSCOPE & OPERATING OTOSCOPE

OPT700.020.030P



PRODUCT DESCRIPTION

Veterinary Diagnostic Set

Introducing a new standard in LED Veterinary Diagnostic Equipment. The ultimate open head design for easier instrument access is a key feature.

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

Set includes:

- Veterinary Operating otoscope
- L28 Practice ophthalmoscope with cobalt blue filter, light sensitive filter and 5 apertures
- 4mm, 5mm and 7mm plastic specula
- 1 x battery handle
- Zip case

Key Features

- TrueTone LED
- Can be 'locked' on to handle
- Excellent LED illumination
- Bayonet tip connection
- Swivelling 3x Magnification lens for enhanced viewing and better instrumentation
- 28 lens from -25d to +40d. (-10 to +10 in single dioptre steps)
- Ultimate Open Head Design Ideal for Instrumentation
- Durable, solid metal construction
- Opticlar and Welch Allyn connection
- Protected fibre optic bundle
- Compatibility with Heine and keeler connectors upon request.
- 6 apertures; large spot, medium spot, small spot, graticule with fixation ring, slit, cobalt blue
- 3 filters; Red free, Clear, Polarised light sensitive
- Dust proofed optics
- Rubber brow rest
- 10 year LED warranty
- 20,000 hour plus LED life guarantees unlimited use with no blown bulb downtime
- 5 year product warranty
- Supplied with 3 plastic tips

Adapt Desk Rechargeable Handle

- Requires 2x c cell batteries
- Robust metal sleeve
- Knurled surface for improved grip
- Rotary control allows the LED brightness to be easily adjusted

Specifications:

- Running duration: > 60 hours of continuous use

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.020.030P
Brand	Opticlar
Diagnostic Head Type	Co-axial Ophthalmoscope, Slit Oscope

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights