

## OPTICLAR DESK DIAGNOSTIC SET WITH OPHTHALMOSCOPE AND VARIOSCOPE

OPT700.020.143



### PRODUCT DESCRIPTION

---

## Veterinary Diagnostic Set - Twin Port Desk Charger

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

#### Set includes:

- Veterinary VarioScope otoscope
- L28 Practice ophthalmoscope with cobalt blue filter, light sensitive filter and 5 apertures
- 2 Adapt rechargeable handle
- Twin port charger
- Set of 3 plastic tips (1 each) of 61.0 x 4.0mm, 61.0 x 6.0mm, 61.0 x 9.0mm
- Zip Case

#### Key Features

- TrueTone LED
- Can be 'locked' on to handle
- Excellent LED illumination
- Bayonet tip connection

- 2x viewing window swivels out of way for easier instrumentation
- 28 lens from -20d to +40d
- Durable, solid metal construction
- Opticlar and Welch Allyn connection
- Compatibility with Heine veterinary tips (Heine Beta 200 Veterinary style)
- 5 apertures: Large, medium and small spot; Graticule and Slit
- 10 year LED warranty
- 20,000 hour plus LED life guarantees unlimited use with no blown bulb downtime
- 5 year product warranty
- Supplied with 3 plastic tips

### Adapt Desk Rechargeable Handle

- Lithium Ion
- Robust metal sleeve
- Knurled surface for maximum grip
- Rotary control allows the LED brightness to be easily adjusted

### Recharging Pod

- Comply with national and international standards
- CE marked
- Single port recharging pod
- Supplied as standard with 3 pin UK plug, other configurations available upon request
- Generous 3m/ 10ft lead

### Specifications:

- Charging time: 2 hours
- Running duration: > 20 hours

### Varioscope compatible with the following plastic specula

- OPT700.000.015
- OPT700.000.016
- OPT700.000.017

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	OPT700.020.143
<b>Brand</b>	Opticlar
<b>Catagory Sub Heading</b>	Save 16%
<b>Diagnostic Head Type</b>	Slit Otoscope

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights