

## PIN INDEX PIPELINE BLOCK - O<sup>2</sup> AND N<sup>2</sup>O

102-85\*BUNDLE



### PRODUCT DESCRIPTION

## Pin Index Pipeline Block

Convert your anaesthetic machine to piped gases using the machines cylinder yokes with the addition of these pin index pipeline blocks.

- For use with ¼ BSP (102-777) and 3/8 (102-778 ) pipeline leads (available separately)
- O<sup>2</sup> and N<sup>2</sup>O version available

### ADDITIONAL INFORMATION

**Animal**

Mixed

**SKU**

102-85\*bundle

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

