

**OPTICLAR METAL SPECULA - MATT BLACK - 57MM X 4.0MM**

OPT700.000.001MB

**PRODUCT DESCRIPTION****Veterinary Metal Specula - Matt Black - 57mm x 4.0mm**

The bayonet fitting is compatible with Opticlar and Heine Slit Otoscope heads. The Matt Black finish completely eliminates glare. Fully autoclaveable.

**Includes:**

- 1 closed metal tip

**Specifications:**

- 57.0mm length
- 4.0mm diameter
- Metal construction
- Autoclaveable
- Bite-resistant

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	OPT700.000.001MB
<b>Brand</b>	Opticlar

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights