

OPTICLAR VISIONVET CLOSED METAL TIP - 150.0MM X 10.0MM

OPT700.000.004

**PRODUCT DESCRIPTION****Veterinary VisionVet Closed Metal Tip Specula - 150mm x 10.0mm diameter**

The bayonet fitting is compatible with Opticlar and Heine Slit Otoloscope heads. Fully autoclaveable

Includes:

- 1 closed metal tip

Specifications:

- 150.0mm length
- 10.0mm diameter

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.000.004
Brand	Opticlar

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights