

OPTICLAR FEMALE CATHETERISATION SPECULUM - 85MM X 9.0MM

OPT700.000.007-2

**PRODUCT DESCRIPTION****Veterinary Female Catheterisation Speculum - 85mm x 9.0mm****For use with Veterinary Slit & Pocket Oscopes**

Designed to aid catheterisation with an overall length of 95.0mm x 9.0mm diameter. The tip has a ventral cut out of approximately 25mm which allows the catheter to be introduced into the urethral orifice under magnification and illumination. Once the catheter has been placed the speculum can be removed and slid back off the catheter (depending upon the size of the catheter) or left on the catheter with the patient. The speculum finish allows easy cleaning whilst the interior has a matt finish to minimise glare and to maximise the amount of light 'thrown' from the tip end. The tip provides excellent visualisation.

Key Features

- Female catheterisation speculum
- Metal construction
- Size 85.0 x 9.0mm catheterisation tip with cutout
- Interior matt finish
- Suitable for both Opticlar and Heine Slit otoscope heads

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.000.007-2
Brand	Opticlar

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights