

OPTICLAR VISIONVET PLASTIC SPECULA - OPEN GUTTER, CLOSED TIP - 90MM X 7.0MM

OPT700.000.025

**PRODUCT DESCRIPTION****Open Gutter, Closed Tips - 90.0mm x 7.0mm****For Veterinary Slit & Pocket Oscopes**

Opticlar VisionVet Closed Plastic Tips

Specifications:

- 90.0mm length
- 7.0mm diameter

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.000.025
Brand	Opticlar

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights